

## Anti-Ribosomal Protein S27 RPS27 Monoclonal Antibody

Catalog Number: M06699S27

### Overview

Product Name	Anti-Ribosomal Protein S27 RPS27 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Ribosomal Protein S27 RPS27 Monoclonal Antibody catalog # M06699S27. Tested in ELISA, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Monoclonal 4A12
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P42677

### Technical Details

Immunogen	Purified recombinant fragment of Ribosomal Protein S27 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ribosomal Protein S27 RPS27 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.