

Anti-Acetyl-Histone H3 (Lys9) H3F3A Rabbit Monoclonal Antibody, Clone#RM161

Catalog Number: M06819-11


Overview

Product Name	Anti-Acetyl-Histone H3 (Lys9) H3F3A Rabbit Monoclonal Antibody, Clone#RM161
Reactive Species	Human, Vertebrates
Description	Boster Bio Anti-Acetyl-Histone H3 (Lys9) H3F3A Rabbit Monoclonal Antibody, Clone#RM161 (Catalog # M06819-11). Tested in WB, ICC, ChIP, IHC, ELISA, Luminex applications. This antibody reacts with Human, Vertebrates.
Application	ChIP, ELISA, IHC, ICC, WB, Luminex
Clonality	Monoclonal RM161
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P84243


Technical Details

Immunogen	An acetyl-peptide corresponding to the Acetyl-Histone H3 (Lys9)
Cross Reactivity	This antibody reacts to Histone H3 acetylated at Lysine 9 (K9ac). No cross-reactivity with other acetylated Lysines in histone H3.
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/mL
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	WB: 0.25 ug/mL - 1 ug/mL ICC: 0.5 ug/mL - 2 ug/mL ChIP: 2 ug/mL - 10 ug/mL IHC: 0.1 ug/mL - 1 ug/mL ELISA: 0.2 ug/mL - 1 ug/mL Luminex: 0.05 ug/mL - 0.5 ug/mL.


Anti-Acetyl-Histone H3 (Lys9) H3F3A Rabbit Monoclonal Antibody, Clone#RM161 (M06819-11) Images




Western Blotting result Western Blot of acid extracts from HeLa cells untreated (-) or treated (+) with sodium butyrate, using RM161 at 0.25 ug/mL, showed a band of histone H3 acetylated at Lysine 9 in treated HeLa.



Specificity Test result RM161 specifically reacts to Histone H3 acetylated at Lysine 9 (K9ac). No cross reactivity with acetylated Lysine 4 (K4ac), Lysine 14 (K14ac), Lysine 18 (K18ac), Lysine 23 (K23ac), Lysine 27 (K27ac), Lysine 36 (K36ac), Lysine 56 (K56ac), lysine 79 (K79ac), or Lysine 122 (K122) in histone H3.





IHC result Immunohistochemistry staining of HepG2 cell using anti-Acetyl-Histone H3 (Lys9) antibody, RM161.



ICC result Immunocytochemistry of HeLa cells treated with sodium butyrate, using Acetyl-Histone H3 (Lys9) Rabbit mAb RM161 (red). Actin filaments have been labeled with fluorescein phalloidin (green).

ELISA result showing specificity Sandwich ELISA against acetylated histone H3 at Lys 9 using HeLa whole cell lysate, treated or untreated with Sodium Butyrate. Using anti-H3CT (RM188, 1 ug/mL) as the capture antibody and biotinylated anti-H3K9ac (RM161, 1 ug/mL) as the detection antibody.



ChIP resultChIP performed on HeLa cells using H3K9ac antibody (RM161, 5ug). Real-time PCR was performed using primers specific to the gene indicated.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Acetyl-Histone H3 (Lys9) H3F3A Rabbit Monoclonal Antibody, Clone#RM161

For Research Use Only. Not for use in diagnostic procedures.