

Anti-Histone H3.3 H3F3A Rabbit Monoclonal Antibody, Clone#RM190

Catalog Number: M06819-29


Overview

Product Name	Anti-Histone H3.3 H3F3A Rabbit Monoclonal Antibody, Clone#RM190
Reactive Species	Vertebrates
Description	Boster Bio Anti-Histone H3.3 H3F3A Rabbit Monoclonal Antibody, Clone#RM190 (Catalog # M06819-29). Tested in WB, ICC, ELISA, Luminex, IF applications. This antibody reacts with Vertebrates.
Application	ELISA, IF, ICC, WB, Luminex
Clonality	Monoclonal RM190
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P84243


Technical Details

Immunogen	A peptide corresponding to human Histone H3.3
Cross Reactivity	This antibody reacts to Histone H3.3, independent of post-translational modifications. No cross-reactivity with Histone H3.1 or other histone proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/mL
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	WB: 0.5 ug/mL - 2 ug/mL ELISA: 0.2 ug/mL - 1 ug/mL Luminex: 0.05 ug/mL - 0.5 ug/mL.

Anti-Histone H3.3 H3F3A Rabbit Monoclonal Antibody, Clone#RM190 (M06819-29) Images



Western Blotting result Western Blot of HeLa whole cell lysate, recombinant Histone H3.3 and Histone H3.1 proteins, using RM190 at 1 ug/mL, showed a band of histone H3.3 in HeLa cells and recombinant H3.3 protein, and no cross activity with Histone H3.1.



Specificity Test result RM190 reacts specifically to Histone H3.3. No cross reactivity with Histone H3.1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H3.3 H3F3A Rabbit Monoclonal Antibody, Clone#RM190

For Research Use Only. Not for use in diagnostic procedures.