

Anti-Trimethyl-Histone H3 (Lys36) H3F3A Rabbit Monoclonal Antibody, Clone#RM155

Catalog Number: M06819-7

Overview

Product Name	Anti-Trimethyl-Histone H3 (Lys36) H3F3A Rabbit Monoclonal Antibody, Clone#RM155
Reactive Species	Human, Vertebrates
Description	Boster Bio Anti-Trimethyl-Histone H3 (Lys36) H3F3A Rabbit Monoclonal Antibody, Clone#RM155 (Catalog # M06819-7). Tested in WB, ChIP, ELISA, Luminex applications. This antibody reacts with Human, Vertebrates.
Application	ChIP, ELISA, WB, Luminex
Clonality	Monoclonal RM155
Formulation	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P84243

Technical Details

Immunogen	A trimethyl-peptide corresponding to trimethyl-Histone H3 (Lys36).
Cross Reactivity	This antibody reacts to Histone H3 trimethylated at Lysine 36 (K36me3). No cross-reactivity with monomethylated Lysine 36 (K36me1) or dimethylated Lysine 36 (K36me2), or other methylations in histone H3.
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/mL
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB: 0.5 ug/mL - 2 ug/mL</p> <p>ChIP: 2 ug/mL- 10 ug/mL</p> <p>ELISA: 0.2 ug/mL - 1 ug/mL</p>

Luminex: 0.1 ug/mL – 0.5 ug/mL.

Anti-Trimethyl-Histone H3 (Lys36) H3F3A Rabbit Monoclonal Antibody, Clone#RM155 (M06819-7) Images




Figure 1. Western Blotting result
Western Blot of recombinant histone H3.3 (1) and acid extracts of HeLa cells (2), using RM155 at 2 ug/mL, showed a band of histone H3 trimethylated at Lysine 36 (K36me3) in HeLa cells.




Figure 2. Specificity Test result
RM155 specifically reacts to Histone H3 trimethylated at Lysine 36 (K36me3). No cross reactivity with non-modified Lysine 36 (K36 Ctrl), monomethylated Lysine 36 (K36me1) or dimethylated Lysine 36 (K36me2), or other methylations in Histone H3.




Figure 3. ChIP result
ChIP performed on HeLa cells using H3K36me3 antibody (RM155, 5ug). Real-time PCR was performed using primers specific to the gene indicated.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Trimethyl-Histone H3 (Lys36) H3F3A Rabbit Monoclonal Antibody, Clone#RM155