

Anti-NEFM/Nf M Rabbit Monoclonal Antibody

Catalog Number: M06821


Overview

| | |
|----------------------|--|
| Product Name | Anti-NEFM/Nf M Rabbit Monoclonal Antibody |
| Reactive Species | Human, Pig, Rat |
| Description | Boster Bio Anti-NEFM/Nf M Rabbit Monoclonal Antibody catalog # M06821. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human, Rat, Pig. |
| Application | IP, IF, IHC, ICC, WB |
| Clonality | Monoclonal AAGI-14 |
| Formulation | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P07197 |


Technical Details

| | |
|---------------------|---|
| Immunogen | A synthesized peptide derived from human NEFM |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5mg/ml |
| Purification | Affinity-chromatography |
| Suggested Dilutions | WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 IP 1:20 |


Anti-NEFM/Nf M Rabbit Monoclonal Antibody (M06821) Images



Western blot analysis of NF-M/NEFM using anti-NF-M/NEFM antibody (M06821). Electrophoresis was performed on a 8% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human 293T whole cell lysates, Lane 2: human SH-SY5Y whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-NF-M/NEFM antigen affinity purified monoclonal antibody (M06821) at 1:500 overnight at 4°C, then washed with TBS-0.1% Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for NF-M/NEFM at approximately 160 kDa. The expected band size for NF-M/NEFM is at 102 kDa.



IHC analysis of NF-M/NEFM using anti-NF-M/NEFM antibody (M06821). NF-M/NEFM was detected in a paraffin-embedded section of rat cerebellum tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-NF-M/NEFM Antibody (M06821) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.



IHC analysis of NF-M/NEFM using anti-NF-M/NEFM antibody (M06821). NF-M/NEFM was detected in a paraffin-embedded section of pig brain tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-NF-M/NEFM Antibody (M06821) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NEFM/Nf M Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.