

Anti-beta II tubulin TUBB2A Monoclonal Antibody

Catalog Number: M06868

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-beta II tubulin TUBB2A Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-beta II tubulin TUBB2A Monoclonal Antibody catalog # M06868. Tested in IP, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, IHC, WB
Clonality	Monoclonal Mix
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q13885

Technical Details


Immunogen	Synthetic Peptide of beta II tubulin
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Suggested Dilutions	WB 1:50000-100000

	IP 1:200 IHC 1:50-300 IF 1:50-200
--	---


Anti-beta II tubulin TUBB2A Monoclonal Antibody (M06868) Images




Western blot analysis of 1) Mouse Brain tissue, 2) Rat Brain tissue, diluted at 1:100000.




1) Input: Mouse Brain Tissue Lysate 2) IP product: IP dilute 1:200 Western blot analysis: primary antibody : TDY191 1:10000 Secondary antibody: Goat anti-Mouse IgG, Light chain specific (S003), 1: 5000



Immunohistochemical analysis of paraffin-embedded Human-colon tissue. 1, beta II tubulin Monoclonal Antibody (Mix) was diluted at 1:200 (4°C, overnight). 2, Sodium citrate pH 6.0 was used for antibody retrieval (>98°C, 20min). 3, Secondary antibody was diluted at 1:200 (room temperature, 30min). Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Rat-kidney tissue. 1, beta II tubulin Monoclonal Antibody (Mix) was diluted at 1:200 (4°C, overnight). 2, Sodium citrate pH 6.0 was used for antibody retrieval (>98°C, 20min). 3, Secondary antibody was diluted at 1:200 (room temperature, 30min). Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Mouse-testis tissue. 1, beta II tubulin Monoclonal Antibody (Mix) was diluted at 1:200 (4°C, overnight). 2, Sodium citrate pH 6.0 was used for antibody retrieval (>98°C, 20min). 3, Secondary antibody was diluted at 1:200 (room temperature, 30min). Negative control was used by secondary antibody only.

Anti-beta II tubulin TUBB2A Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.