

Anti-TC10 Rabbit Monoclonal Antibody

Catalog Number: M06900-1

About RHOQ

Putative transcription factor involved in pancreas development and function.


Overview

Product Name	Anti-TC10 Rabbit Monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-TC10 Rabbit Monoclonal Antibody catalog # M06900-1. Tested in WB, ICC/IF applications. This antibody reacts with Human, Rat.
Application	IF, ICC, WB
Clonality	Monoclonal 21R60
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P17081

Technical Details

Immunogen	A synthesized peptide derived from human TC10
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB 1:500-1:2000 ICC/IF 1:50-1:200

Anti-TC10 Rabbit Monoclonal Antibody (M06900-1) Images



Western blot analysis of TC10 expression in 293T cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TC10 Rabbit Monoclonal Antibody