

Anti-RNF114 Antibody (Monoclonal, 34R78)

Catalog Number: M06969

About RNF114

Predicted to enable ubiquitin protein ligase activity. Predicted to be involved in protein polyubiquitination and ubiquitin-dependent protein catabolic process. Located in cytosol and plasma membrane. Biomarker of male infertility.

Overview

Product Name	Anti-RNF114 Antibody (Monoclonal, 34R78)
Reactive Species	Human
Description	Boster Bio Anti-RNF114 Antibody (Monoclonal, 34R78) catalog # M06969. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 34R78
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y508

Technical Details

Immunogen	Recombinant protein within human RNF114 aa 3-226.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RNF114 Antibody (Monoclonal, 34R78)

For Research Use Only. Not for use in diagnostic procedures.