

## Anti-SWAP70 monoclonal Antibody

Catalog Number: M06972


### Overview

Product Name	Anti-SWAP70 monoclonal Antibody
Reactive Species	Hamster, Human, Mouse, Rat
Description	Boster Bio Anti-SWAP70 monoclonal Antibody catalog # M06972. Tested in WB,IP applications. This antibody reacts with Human,Mouse,Rat,Hamster.
Application	IP, WB
Clonality	Monoclonal 1D12
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UH65


### Technical Details

Immunogen	Recombinant full length Human SWAP70.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1000 IP: 1:50-100

## Anti-SWAP70 monoclonal Antibody (M06972) Images



Western blot detection of SWAP70 in HeLa, 3T3, CHO-K1 and C6 cell lysates and using SWAP70 antibody (1:500 diluted) .



Immunoprecipitation analysis of HeLa cell lysate using SWAP70 antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SWAP70 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.