

Anti-TBL1 (4H2) monoclonal TBL1X Antibody

Catalog Number: M07014


Overview

Product Name	Anti-TBL1 (4H2) monoclonal TBL1X Antibody
Reactive Species	Human
Description	Boster Bio Anti-TBL1 (4H2) monoclonal TBL1X Antibody catalog # M07014. Tested in WB,ICC applications. This antibody reacts with Human.
Application	ICC, WB
Clonality	Monoclonal 1D12
Formulation	Mouse IgG2b, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O60907


Technical Details

Immunogen	Purified recombinant human TBL1 protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000 ICC: 1:100-500

Anti-TBL1 (4H2) monoclonal TBL1X Antibody (M07014) Images



Western blot detection of TBL1X in SW480 cell lysates using TBL1X mouse mAb (1:1000 diluted).



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-TBL1X mouse mAb (dilution 1:100).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBL1 (4H2) monoclonal TBL1X Antibody

For Research Use Only. Not for use in diagnostic procedures.