

Anti-Erlin-2/ERLIN2 Antibody Picoband® (monoclonal, 3H9A2) PE Conjugated

Catalog Number: M07042-1-PE

About ERLIN2

Erlin-2 is a protein that in humans is encoded by the ERLIN2 gene. This gene encodes a member of the SPFH domain-containing family of lipid raft-associated proteins. The encoded protein is localized to lipid rafts of the endoplasmic reticulum and plays a critical role in inositol 1,4,5-trisphosphate (IP3) signaling by mediating ER-associated degradation of activated IP3 receptors. Mutations in this gene are a cause of spastic paraplegia-18 (SPG18). Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Erlin-2/ERLIN2 Antibody Picoband® (monoclonal, 3H9A2) PE Conjugated |
| Reactive Species | Human |
| Application | Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Monoclonal 3H9A2 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Mouse |
| Uniprot ID | O94905 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human Erlin-2/ERLIN2 recombinant protein (Position: D87-N339). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG2b |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Erlin-2/ERLIN2 Antibody (monoclonal, 3H9A2) - PE

For Research Use Only. Not for use in diagnostic procedures.