

Anti-SMYD1 monoclonal Antibody

Catalog Number: M07134


Overview

Product Name	Anti-SMYD1 monoclonal Antibody
Reactive Species	Rat
Description	Boster Bio Anti-SMYD1 monoclonal Antibody catalog # M07134. Tested in WB applications. This antibody reacts with Rat.
Application	WB
Clonality	Monoclonal Clone: 9:4-7 IgG2b kappa
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q8NB12

Technical Details

Immunogen	Recombinant full length Human SMYD1.
Isotype	IgG2b kappa
Form	Liquid
Concentration	1.0 mg/mL by UV absorbance at 280 nm
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-2000

Anti-SMYD1 monoclonal Antibody (M07134) Images



Western blot detection of SMYD1 in rat skeletal muscle lysates and using SMYD1 antibody (1:2000 diluted) .

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SMYD1 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.