

Anti-MAGEA4 Mouse Monoclonal Antibody [Clone ID: OTI2C1]

Catalog Number: M07175


Overview

Product Name	Anti-MAGEA4 Mouse Monoclonal Antibody [Clone ID: OTI2C1]
Reactive Species	Human
Description	Boster Bio MAGEA4 mouse monoclonal antibody, clone OTI2C1 (formerly 2C1). Catalog# M07175. Tested in IF, IHC, WB. This antibody reacts with Human.
Application	IF, IHC, WB
Clonality	Monoclonal OTI2C1
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P43358


Technical Details

Immunogen	Full length human recombinant protein of human MAGEA4 (NP_001011550) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)


Anti-MAGEA4 Mouse Monoclonal Antibody [Clone ID: OTI2C1] (M07175) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MAGEA4 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MAGEA4.




Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-MAGEA4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-MAGEA4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-MAGEA4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human colon tissue within the normal limits using anti-MAGEA4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Anti-MAGEA4 mouse monoclonal antibody (M07175)
immunofluorescent staining of COS7 cells transiently
transfected by pCMV6-ENTRY MAGEA4.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MAGEA4 Mouse Monoclonal Antibody [Clone ID: OTI2C1]