

Anti-SULT1C2 Mouse Monoclonal Antibody [Clone ID: OTI5H1]

Catalog Number: M07212-1


Overview

Product Name	Anti-SULT1C2 Mouse Monoclonal Antibody [Clone ID: OTI5H1]
Reactive Species	Dog, Human, Monkey, Mouse
Description	Boster Bio SULT1C2 mouse monoclonal antibody, clone OTI5H1 (formerly 5H1). Catalog# M07212-1. Tested in WB. This antibody reacts with Human, Monkey, Mouse, Dog.
Application	WB
Clonality	Monoclonal OTI5H1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O00338


Technical Details

Immunogen	Full length human recombinant protein of human SULT1C2 (NP_001047) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1.12 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500~2000

Anti-SULT1C2 Mouse Monoclonal Antibody [Clone ID: OTI5H1] (M07212-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SULT1C2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SULT1C2.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-SULT1C2 monoclonal antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SULT1C2 Mouse Monoclonal Antibody [Clone ID: OTI5H1]

For Research Use Only. Not for use in diagnostic procedures.