

Anti-IL18BP monoclonal Antibody

Catalog Number: M07261


Overview

Product Name	Anti-IL18BP monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL18BP monoclonal Antibody catalog # M07261. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 3C10F7 / 3C10G5 / PMC01 / 4D12B3 / 3C10F7C11
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O95998

Technical Details

Immunogen	Recombinant full length Human IL18BP.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:400-1000

Anti-IL18BP monoclonal Antibody (M07261) Images



Western blot detection of IL18BP antibody in CHO-K1 cells transfected by IL18BP-PDGFR fusion protein and K562 cell lysates using IL18BP antibody (1:400-1:1000 diluted) .

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL18BP monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.