

Anti-Histone H2B (mono methyl R79) HIST2H2BE Rabbit Monoclonal Antibody

Catalog Number: M07286

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Histone H2B (mono methyl R79) HIST2H2BE Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Histone H2B (mono methyl R79) HIST2H2BE Rabbit Monoclonal Antibody catalog # M07286. Tested in WB application. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal DFD-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16778

Technical Details

Immunogen	A synthesized peptide derived from human Histone H2B (mono methyl R79)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml

Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

Anti-Histone H2B (mono methyl R79) HIST2H2BE Rabbit Monoclonal Antibody (M07286) Images



Western blot analysis of Histone H2B (mono methyl R79) expression in HeLa cell lysate.

[Anti-Histone H2B \(mono methyl R79\) HIST2H2BE Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.