

## Anti-Fibulin-2 FBLN2 Monoclonal Antibody

Catalog Number: M07304


### Overview

Product Name	Anti-Fibulin-2 FBLN2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Fibulin-2 FBLN2 Monoclonal Antibody catalog # M07304. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 7H4F6/3B9A11
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P98095

### Technical Details

Immunogen	Purified recombinant fragment of Fibulin-2 (aa180-440) expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-Fibulin-2 FBLN2 Monoclonal Antibody (M07304) Images



Western Blot analysis using Fibulin-2 Monoclonal Antibody against truncated FBLN2-Trx recombinant protein (1) and truncated FBLN2 (aa28-444)-hlgGfc transfected COS7 cell lysate (2).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fibulin-2 FBLN2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.