

Anti-COG1 Antibody (C-term)

Catalog Number: M07344

About COG1

Required for normal Golgi function (By similarity).


Overview

Product Name	Anti-COG1 Antibody (C-term)
Reactive Species	Human
Description	Boster Bio Anti-COG1 Antibody (C-term) (Catalog # M07344). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WTW3

Technical Details

Immunogen	This COG1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 887-920 amino acids from the C-terminal region of human COG1.
Predicted Reactive Species	Mouse
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-COG1 Antibody (C-term) (M07344) Images



Western blot analysis of lysate from A549 cell line, using COG1 Antibody (C-term). M07344 was diluted at 1:1000. A goat anti-rabbit IgG H&L (HRP) at 1:10000 dilution was used as the secondary antibody. Lysate at 20ug.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COG1 Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.