

## Anti-HNF6 (ONECUT1) Mouse Monoclonal Antibody [Clone ID: OTI4F12]

Catalog Number: M07349


### Overview

Product Name	Anti-HNF6 (ONECUT1) Mouse Monoclonal Antibody [Clone ID: OTI4F12]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio HNF6 mouse monoclonal antibody, clone OTI4F12 (formerly 4F12). Catalog# M07349. Tested in IF, IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Monoclonal OTI4F12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9UBC0


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 85-198 of human HNF6 (NP_004489) produced in E.coli.
Isotype	IgG1
Concentration	0.72 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IHC 1:50 IF 1:100


## Anti-HNF6 (ONECUT1) Mouse Monoclonal Antibody [Clone ID: OTI4F12] (M07349) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HNF6 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HNF6.



Immunohistochemical staining of paraffin-embedded Carcinoma of Human liver tissue using anti-HNF6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Anti-HNF6 mouse monoclonal antibody (M07349) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HNF6.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HNF6 (ONECUT1) Mouse Monoclonal Antibody [Clone ID: OTI4F12]

For Research Use Only. Not for use in diagnostic procedures.