

Anti-Matrilin-1 MATN1 Monoclonal Antibody

Catalog Number: M07385


Overview

Product Name	Anti-Matrilin-1 MATN1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Matrilin-1 MATN1 Monoclonal Antibody catalog # M07385. Tested in ELISA, IF, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 5A8
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P21941


Technical Details

Immunogen	Purified recombinant fragment of human Matrilin-1 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-Matrilin-1 MATN1 Monoclonal Antibody (M07385) Images



Western Blot analysis using Matrilin-1 Monoclonal Antibody against HEK293 (1) and MATN1-hlgGfC transfected HEK293 (2) cell lysate.



Immunofluorescence analysis of HeLa cells using Matrilin-1 Monoclonal Antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Matrilin-1 MATN1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.