

Anti-Akt1 PXX Monoclonal Antibody

Catalog Number: M07555


Overview

Product Name	Anti-Akt1 PXX Monoclonal Antibody
Reactive Species	Mouse
Description	Boster Bio Anti-Akt1 PXX Monoclonal Antibody catalog # M07555. Tested in ELISA, WB applications. This antibody reacts with Mouse.
Application	WB
Clonality	Monoclonal 2D9-D2-H3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q7Z7A4


Technical Details

Immunogen	Purified recombinant human AKT1 protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:2000

Anti-Akt1 PXX Monoclonal Antibody (M07555) Images



Western blotting validation for Anti-Akt1 PXX Monoclonal Antibody M07555 Western blot detection of AKT1 in 3T3 cell lysates using AKT1 mouse mAb (dilution 1:2000)



Western blotting validation for Anti-Akt1 PXX Monoclonal Antibody M07555 Western blot detection of AKT1 in CHO-K1 cell lysate (B) and CHO-K1 transfected by AKT1 (A) cell lysate using AKT1 mouse mAb (1:1000 diluted). Predicted band size: 60kDa. Observed band size: 60kDa. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Akt1 PXX Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.