

Anti-RFX3 Mouse Monoclonal Antibody [Clone ID: OTI3F7]

Catalog Number: M07674


Overview

Product Name	Anti-RFX3 Mouse Monoclonal Antibody [Clone ID: OTI3F7]
Reactive Species	Human, Mouse
Description	Boster Bio RFX3 mouse monoclonal antibody, clone OTI3F7. Catalog# M07674. Tested in WB. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal OTI3F7
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P48380

Technical Details

Immunogen	Full length human recombinant protein of human RFX3 (NP_602304) produced in HEK293T cell.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/

Anti-RFX3 Mouse Monoclonal Antibody [Clone ID: OTI3F7] (M07674) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RFX3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RFX3.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RFX3 Mouse Monoclonal Antibody [Clone ID: OT13F7]