

## Anti-ZNF365 Mouse Monoclonal Antibody [Clone ID: OTI8G7]

Catalog Number: M07837


### Overview

Product Name	Anti-ZNF365 Mouse Monoclonal Antibody [Clone ID: OTI8G7]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio ZNF365 mouse monoclonal antibody, clone OTI8G7 (formerly 8G7). Catalog# M07837. Tested in FC, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal OTI8G7
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q70YC4


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 108-407 of human ZNF365 (NP_055766) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 Flow Cytometry 1:100

## Anti-ZNF365 Mouse Monoclonal Antibody [Clone ID: OTI8G7] (M07837) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF365 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF365.



HEK293T cells transfected with either ZNF365 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ZNF365 antibody (M07837)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZNF365 Mouse Monoclonal Antibody [Clone ID: OTI8G7]

For Research Use Only. Not for use in diagnostic procedures.