

## Anti-TFEC Mouse Monoclonal Antibody [Clone ID: OTI1B4]

Catalog Number: M07841


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-TFEC Mouse Monoclonal Antibody [Clone ID: OTI1B4]   |
| Reactive Species     | Human  |
| Description          | Boster Bio TFEC mouse monoclonal antibody, clone OTI1B4. Catalog# M07841. Tested in WB. This antibody reacts with Human.   |
| Application          | WB   |
| Clonality            | Monoclonal OTI1B4  |
| Formulation          | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received.  |
| Host                 | Mouse  |
| Uniprot ID           | O14948   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Full length human recombinant protein of human TFEC (NP_001018068) produced in E.coli.                    |
| Isotype             | IgG1  |
| Concentration       | 1 mg/ml   |
| Purification        | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:2000   |

## Anti-TFEC Mouse Monoclonal Antibody [Clone ID: OTI1B4] (M07841) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TFEC (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TFEC (1:2000).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TFEC Mouse Monoclonal Antibody [Clone ID: OTI1B4]

For Research Use Only. Not for use in diagnostic procedures.