

Anti-SLC44A4 Reference Antibody (ASG-5ME)

Catalog Number: M08008

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-SLC44A4 Reference Antibody (ASG-5ME)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SLC44A4 Reference Antibody (ASG-5ME) (Catalog # M08008). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.936EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q53GD3


Technical Details

Isotype	IgG2SA
Form	Liquid

Anti-SLC44A4 Reference Antibody (ASG-5ME) (M08008) Images



Anti-SLC44A4 Reference Antibody (ASG-5ME) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-SLC44A4 Reference Antibody (ASG-5ME) is more than 95%

Anti-SLC44A4 Reference Antibody (ASG-5ME)

For Research Use Only. Not for use in diagnostic procedures.