

Anti-CHRDL1 monoclonal Antibody

Catalog Number: M08014

Overview

Product Name	Anti-CHRDL1 monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CHRDL1 monoclonal Antibody catalog # M08014. Tested in WB,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	ICC, WB
Clonality	Monoclonal Clone: 22C5.H7.A6 IgG2a kappa
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9BU40

Technical Details

Immunogen	Recombinant full length Human CHRDL1.
Isotype	IgG2a kappa
Form	Liquid
Concentration	1.0 mg/mL by UV absorbance at 280 nm
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000 ICC: 1:300

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CHRDL1 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.