

Anti-Mast Cell Tryptase Rabbit Monoclonal Antibody

Catalog Number: M08039


Overview

| | |
|----------------------|--|
| Product Name | Anti-Mast Cell Tryptase Rabbit Monoclonal Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Mast Cell Tryptase Rabbit Monoclonal Antibody catalog # M08039. Tested in WB, IHC, IP applications. This antibody reacts with Human, Mouse, Rat. |
| Application | IP, IHC, WB |
| Clonality | Monoclonal FHH-20 |
| Formulation | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q15661 |


Technical Details

| | |
|---------------------|--|
| Immunogen | A synthesized peptide derived from human Mast Cell Tryptase |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Purification | Affinity-chromatography |
| Suggested Dilutions | WB 1:500-1:2000 IHC 1:50-1:200 IP 1:50 |

Anti-Mast Cell Tryptase Rabbit Monoclonal Antibody (M08039) Images



Western blot analysis of Mast Cell Tryptase expression in Human tonsil lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using Mast Cell Tryptase Antibody.

1 Publications Citing This Product

1. PubMed ID: 30090343, Acute exposure of ozone induced pulmonary injury and the protective role of vitamin E through the Nrf2 pathway in Balb/c mice

Visit bosterbio.com/anti-mast-cell-tryptase-rabbit-monoclonal-antibody-m08039-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mast Cell Tryptase Rabbit Monoclonal Antibody