

Anti-Rex1 (ZFP42) Mouse Monoclonal Antibody [Clone ID: OTI3H9]

Catalog Number: M08177


Overview

Product Name	Anti-Rex1 (ZFP42) Mouse Monoclonal Antibody [Clone ID: OTI3H9]
Reactive Species	Human
Description	Boster Bio ZFP42 mouse monoclonal antibody, clone OTI3H9 (formerly 3H9). Catalog# M08177. Tested in WB. This antibody reacts with Human.
Conjugate	Unconjugated
Application	WB
Clonality	Monoclonal OTI3H9
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q96MM3

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1-235 of human ZFP42 (NP_777560) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:2000

Anti-Rex1 (ZFP42) Mouse Monoclonal Antibody [Clone ID: OTI3H9] (M08177) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZFP42 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZFP42 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rex1 (ZFP42) Mouse Monoclonal Antibody [Clone ID: OTI3H9]

For Research Use Only. Not for use in diagnostic procedures.