

## Anti-CD11 Reference Antibody (Rovelizumab)

Catalog Number: M08209


### Overview

Product Name	Anti-CD11 Reference Antibody (Rovelizumab)
Reactive Species	Human
Description	Boster Bio Anti-CD11 Reference Antibody (Rovelizumab) (Catalog # M08209). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.949EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q9UQ88


### Technical Details

Isotype	IgG4
Form	Liquid

## Anti-CD11 Reference Antibody (Rovelizumab) (M08209) Images



Anti-CD11 Reference Antibody (Rovelizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-CD11 Reference Antibody (Rovelizumab) is more than 99.22%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD11 Reference Antibody (Rovelizumab)

For Research Use Only. Not for use in diagnostic procedures.