

Anti-USP24 Rabbit Monoclonal Antibody

Catalog Number: M08241


Overview

| | |
|----------------------|--|
| Product Name | Anti-USP24 Rabbit Monoclonal Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-USP24 Rabbit Monoclonal Antibody catalog # M08241. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Monoclonal 26U83 |
| Formulation | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9UPU5 |

Technical Details

| | |
|---------------------|--|
| Immunogen | A synthesized peptide derived from human USP24 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Purification | Affinity-chromatography |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB 1:500-1:2000</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p> |

Anti-USP24 Rabbit Monoclonal Antibody (M08241) Images



Western blot analysis of USP24 expression in HeLa cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-USP24 Rabbit Monoclonal Antibody