

Anti-Human IgG3 IGHG3 Monoclonal Antibody

Catalog Number: M08374


Overview

Product Name	Anti-Human IgG3 IGHG3 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Human IgG3 IGHG3 Monoclonal Antibody catalog # M08374. Tested in WB, IP applications. This antibody reacts with Human.
Application	IP, WB
Clonality	Monoclonal ACIA-9
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01860

Technical Details

Immunogen	A synthesized peptide derived from human Human IgG3 IgG is a monomeric immunoglobulin, built of two heavy chains gamma and two light chains. Each molecule has two antigen binding sites. This is the most abundant immunoglobulin and is approximately equally distributed in blood and in tissue liquids, constituting 75% of serum immunoglobulins in humans. There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%) .
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IP 1:50

Anti-Human IgG3 IGHG3 Monoclonal Antibody (M08374) Images



Western blot analysis of Human IgG3 expression in human plasma lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Human IgG3 IGHG3 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.