

## Anti-LYPD3 / C4.4A Reference Antibody (Iupartumab)

Catalog Number: M08396

### Overview

Product Name	Anti-LYPD3 / C4.4A Reference Antibody (Iupartumab)
Reactive Species	Cynomolgus monkey, Human, Mouse
Description	Boster Bio Anti-LYPD3 / C4.4A Reference Antibody (Iupartumab) (Catalog # M08396). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.950EU/ug, determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O95274

### Technical Details

Isotype	IgG1
Form	Liquid

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LYPD3 / C4.4A Reference Antibody (Iupartumab)

For Research Use Only. Not for use in diagnostic procedures.