

Anti-LYPD3 Antibody (Monoclonal, 32L82)

Catalog Number: M08396-2

About LYPD3

Ly6/PLAUR domain-containing protein 3 is a protein that in humans is encoded by the LYPD3 gene. Ly6 / PLAUR domain-containing protein 3, also known as GPI-anchored metastasis-associated protein C4.4A homolog, Matrigel-induced gene C4 protein, MIG-C4, and LYPD3, is a cell membrane protein that contains two UPAR/Ly6 domains. Human LYPD3 contains two UPAR/Ly6 domains. LYPD3 is expressed in the placenta, skin, and urothelium. It is found in suprabasal keratinocytes of chronic wounds. Weak expression of LYPD3 is found in the esophagus and peripheral blood mononuclear cells. It is found in the majority of primary and metastatic transitional cell carcinomas (TCCs) and as well in breast cancer tissues, but not in adjacent normal tissues. High expression of LYPD3 is found in the tumor component of some noninvasive superficial lesions and invasive and metastatic urothelial cancers. LYPD3 is up-regulated in migrating keratinocytes during epithelisation of incisional skin wounds. LYPD3 supports cell migration. It may be involved in urothelial cell-matrix interactions. It may also be involved in tumor progression.

Overview

Product Name	Anti-LYPD3 Antibody (Monoclonal, 32L82)
Reactive Species	Human
Description	Boster Bio Anti-LYPD3 Antibody (Monoclonal, 32L82) catalog # M08396-2. Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal 32L82
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95274

Technical Details

Immunogen	Synthetic peptide.
Form	Liquid
Concentration	500 ug/ml

Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LYPD3 Antibody (Monoclonal, 32L82)

For Research Use Only. Not for use in diagnostic procedures.