

Anti-PMF1 Rabbit Monoclonal Antibody

Catalog Number: M08465

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-PMF1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PMF1 Rabbit Monoclonal Antibody catalog # M08465. Tested in WB, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, WB
Clonality	Monoclonal 30P63
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6P1K2

Technical Details

Immunogen	A synthesized peptide derived from human PMF1
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

	IP 1:50
--	---------

Anti-PMF1 Rabbit Monoclonal Antibody (M08465) Images



Western blot analysis of PMF1 expression in (1) HeLa cell lysate; (2) RAW 264.7 cell lysate.

Anti-PMF1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.