

Anti-FUNDC1 Antibody (Monoclonal, 32F72)

Catalog Number: M08688

About FUNDC1

This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known.

Overview

Product Name	Anti-FUNDC1 Antibody (Monoclonal, 32F72)
Reactive Species	Human, Rat
Description	Boster Bio Anti-FUNDC1 Antibody (Monoclonal, 32F72) catalog # M08688. Tested in WB, IHC applications. This antibody reacts with Human, Rat.
Application	IHC, WB
Clonality	Monoclonal 32F72
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8IVP5

Technical Details

Immunogen	Recombinant protein within N terminal human FUNDC1.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FUNDC1 Antibody (Monoclonal, 32F72)

For Research Use Only. Not for use in diagnostic procedures.