

Anti-FUBP3 Antibody (N-term)

Catalog Number: M08695

About FUBP3

May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-FUBP3 Antibody (N-term)
Reactive Species	Human
Description	Boster Bio Anti-FUBP3 Antibody (N-term) (Catalog # M08695). Tested in IF, WB application(s). This antibody reacts with Human.
Application	IF, WB
Clonality	Monoclonal 1216CT820.155.150
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q96I24


Technical Details

Immunogen	Purified His-tagged FUBP3 protein was used to produced this monoclonal antibody.
Predicted Reactive Species	Bovine
Isotype	IgG1
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	IF: 1:25 WB: 1:2000


Anti-FUBP3 Antibody (N-term) (M08695) Images



Fluorescent image of HeLa cells stained with FUBP3 Antibody (N-term). FUBP3 Antibody was diluted at 1:25 dilution. An Alexa Fluor® 488-conjugated goat anti-mouse IgG at 1:400 dilution was used as the secondary antibody (green). Cytoplasmic actin was counterstained with Alexa Fluor® 555 conjugated with Phalloidin (red).



Fluorescent image of HeLa cells stained with FUBP3 Antibody (N-term). FUBP3 Antibody was diluted at 1:25 dilution. An Alexa Fluor® 488-conjugated goat anti-mouse IgG at 1:400 dilution was used as the secondary antibody (green). Cytoplasmic actin was counterstained with Alexa Fluor® 555 conjugated with Phalloidin (red).



FUBP3 Antibody (N-term) western blot analysis in HeLa,U251 cell line lysates (35ug/lane).This demonstrates the FUBP3 antibody detected the FUBP3 protein (arrow).

Anti-FUBP3 Antibody (N-term)

For Research Use Only. Not for use in diagnostic procedures.