

## Anti-TEM7R / PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R)

Catalog Number: M08798


### Overview

Product Name	Anti-TEM7R / PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TEM7R / PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R) (Catalog # M08798). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.951EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q6UX71


### Technical Details

Isotype	IgG1
Form	Liquid

## Anti-TEM7R / PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R) (M08798) Images



Anti-TEM7R/PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TEM7R/PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R) is more than 95%

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TEM7R / PLXDC2 Reference Antibody (Fox Chase patent anti-TEM7R)

For Research Use Only. Not for use in diagnostic procedures.