

Anti-Histone H1.0 H1F0 Rabbit Monoclonal Antibody

Catalog Number: M08821


Overview

| | |
|----------------------|--|
| Product Name | Anti-Histone H1.0 H1F0 Rabbit Monoclonal Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Histone H1.0 H1F0 Rabbit Monoclonal Antibody catalog # M08821. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human, Mouse, Rat. |
| Application | IF, IHC, ICC, WB |
| Clonality | Monoclonal HAG-8 |
| Formulation | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P07305 |


Technical Details

| | |
|---------------------|---|
| Immunogen | A synthesized peptide derived from human Histone H1.0 |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5mg/ml |
| Purification | Affinity-chromatography |
| Suggested Dilutions | WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 |


Anti-Histone H1.0 H1F0 Rabbit Monoclonal Antibody (M08821) Images



Western blot analysis of Histone H1.0 expression in Human kidney lysate.



Immunohistochemical analysis of paraffin-embedded human testis, using Histone H1.0 Antibody.



Immunofluorescent analysis of HepG2 cells, using Histone H1.0 Antibody.

1 Publications Citing This Product

1. PubMed ID: 24503896, Role of endoplasmic reticulum stress in apoptosis of differentiated mouse podocytes induced by high glucose

Visit bosterbio.com/anti-histone-h1-0-rabbit-monoclonal-antibody-m08821-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H1.0 H1F0 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.