

Anti-ASRGL1 (Asparaginase Like 1) Monoclonal Antibody

Catalog Number: M08869

Overview

Product Name	Anti-ASRGL1 (Asparaginase Like 1) Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ASRGL1 (Asparaginase Like 1) Monoclonal Antibody catalog # M08869. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Monoclonal CRASH/1289
Formulation	Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.
Storage Instructions	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Host	Mouse
Uniprot ID	Q7L266

Technical Details

Immunogen	Recombinant human full-length ASRGL1
Isotype	IgG1, kappa
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A/G
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used Immunohistology (Formalin-fixed) (0.5-1.0ug/ml for 30 minutes at RT) Immunofluorescence (0.5-1ug/ml)</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-ASRGL1 (Asparaginase Like 1) Monoclonal Antibody (M08869) Images



Formalin-fixed, paraffin-embedded human Endometrium stained with ASRGL1 Monoclonal Antibody (CRASH/1289).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ASRGL1 (Asparaginase Like 1) Monoclonal Antibody