

Anti-TRIM9 Mouse Monoclonal Antibody [Clone ID: OTI2D6]

Catalog Number: M08967


Overview

Product Name	Anti-TRIM9 Mouse Monoclonal Antibody [Clone ID: OTI2D6]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio TRIM9 mouse monoclonal antibody, clone OTI2D6 (formerly 2D6). Catalog# M08967. Tested in FC, IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI2D6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9C026


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 284-669 of human TRIM9 (NP_443210) produced in E.coli.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IF 1:100 Flow Cytometry 1:100


Anti-TRIM9 Mouse Monoclonal Antibody [Clone ID: OTI2D6] (M08967) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRIM9 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRIM9.



Western blot analysis of extracts (10ug) from 2 different cell lines by using anti-TRIM9 monoclonal antibody at 1:200 dilution.



HEK293T cells transfected with either TRIM9 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-TRIM9 antibody (M08967)



Anti-TRIM9 mouse monoclonal antibody (M08967) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TRIM9.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.