

Anti-Histone H2B Rabbit Monoclonal Antibody, Clone#RM230

Catalog Number: M08989


Overview

Product Name	Anti-Histone H2B Rabbit Monoclonal Antibody, Clone#RM230
Reactive Species	Vertebrates
Description	Boster Bio Anti-Histone H2B Rabbit Monoclonal Antibody, Clone#RM230 (Catalog # M08989). Tested in WB, ICC, ELISA, Luminex applications. This antibody reacts with Vertebrates.
Application	ELISA, ICC, WB, Luminex
Clonality	Monoclonal RM230
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P06899


Technical Details

Immunogen	A peptide corresponding to the C-terminus of human Histone H2B.
Cross Reactivity	This antibody reacts to the Histone H2B protein, independent of post-translational modifications. No cross-reactivity with other histone proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	1 mg/mL
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	WB: 0.1 ug/mL - 1 ug/mL ICC: 0.5 ug/mL - 1 ug/mL ELISA: 0.2 ug/mL - 1 ug/mL Luminex: 0.2 ug/mL - 1 ug/mL.

Anti-Histone H2B Rabbit Monoclonal Antibody, Clone#RM230 (M08989) Images



Western Blotting result Western Blot of recombinant Histone H2A, H2B, the whole cell lysates of HeLa, A375, SK-MEL-2, A431, and K562, using RM230 at 0.2 ug/mL



ICC result Immunocytochemical staining of HeLa cells, using Histone H2B Rabbit mAb RM230 (red). Actin filaments have been labeled with fluorescein phalloidin (green).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone H2B Rabbit Monoclonal Antibody, Clone#RM230

For Research Use Only. Not for use in diagnostic procedures.