

## Anti-LIN7 (LIN7B) Mouse Monoclonal Antibody [Clone ID: OT11C9]

Catalog Number: M09113


### Overview

Product Name	Anti-LIN7 (LIN7B) Mouse Monoclonal Antibody [Clone ID: OT11C9]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-LIN7B mouse monoclonal antibody, clone OT11C9 (formerly 1C9). Catalog# M09113. Tested in FC, IF, IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal OT11C9
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9HAP6


### Technical Details

Immunogen	Full length human recombinant protein of human LIN7B (NP_071448) produced in HEK293T cell.
Isotype	IgG1
Concentration	0.58 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150 IF 1:100 Flow Cytometry 1:100


## Anti-LIN7 (LIN7B) Mouse Monoclonal Antibody [Clone ID: OTI1C9] (M09113) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LIN7B (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LIN7B.




Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-LIN7B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-LIN7B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-LIN7B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



HEK293T cells transfected with either LIN7B (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-LIN7B antibody (M09113)



Anti-LIN7B mouse monoclonal antibody (M09113)  
immunofluorescent staining of COS7 cells transiently  
transfected by pCMV6-ENTRY LIN7B.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LIN7 (LIN7B) Mouse Monoclonal Antibody [Clone ID: OT11C9]

For Research Use Only. Not for use in diagnostic procedures.