

Anti-HAO1 Mouse Monoclonal Antibody [Clone ID: OTI5A1]

Catalog Number: M09159-2


Overview

| | |
|----------------------|--|
| Product Name | Anti-HAO1 Mouse Monoclonal Antibody [Clone ID: OTI5A1] |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio HAO1 mouse monoclonal antibody, clone OTI5A1 (formerly 5A1). Catalog# M09159-2. Tested in FC, IF, WB. This antibody reacts with Human, Mouse, Rat. |
| Application | Flow Cytometry, IF, WB |
| Clonality | Monoclonal OTI5A1 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C as received. |
| Host | Mouse |
| Uniprot ID | Q9UJM8 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Full length human recombinant protein of human HAO1 (NP_060015) produced in HEK293T cell. |
| Isotype | IgG1 |
| Concentration | 0.68 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB: 1:500-1:2000 IF: 1:100 Flow Cytometry: 1:50-1:100 |


Anti-HAO1 Mouse Monoclonal Antibody [Clone ID: OTI5A1] (M09159-2) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HAO1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HAO1.



HEK293T cells transfected with either HAO1 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HAO1 antibody (M09159-2)



Anti-HAO1 mouse monoclonal antibody (M09159-2) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HAO1.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HAO1 Mouse Monoclonal Antibody [Clone ID: OTI5A1]

For Research Use Only. Not for use in diagnostic procedures.