

Anti-TEAD2 Antibody (Monoclonal, 32T77)

Catalog Number: M09168

About TEAD2

TEAD2 (ETF, ETEF-1, TEF-4), together with TEAD1, defines a novel family of transcription factors, the TEAD family, highly conserved through evolution. TEAD proteins were notably found in *Drosophila* (Scalloped), *C. elegans* (egl-44), *S. cerevisiae* and *A. nidulans*. TEAD2 has been less studied than TEAD1 but a few studies revealed its role during development. TEAD2 is a member of the mammalian TEAD transcription factor family (initially named the transcriptional enhancer factor (TEF) family), which contain the TEA/ATTS DNA-binding domain.[7] Members of the family in mammals are TEAD1, TEAD2, TEAD3, TEAD4.

Overview

Product Name	Anti-TEAD2 Antibody (Monoclonal, 32T77)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TEAD2 Antibody (Monoclonal, 32T77) catalog # M09168. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal 32T77
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15562

Technical Details

Immunogen	Recombinant protein within human TEAD2 aa 103-298.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TEAD2 Antibody (Monoclonal, 32T77)

For Research Use Only. Not for use in diagnostic procedures.