

Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab)

Catalog Number: M09198-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab) (Catalog # M09198-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.958EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P55899


Technical Details

Isotype	IgG4SP
Form	Liquid


Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab) (M09198-2) Images



Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab) is more than 97.55%



Immobilized FcRn(FCGRT) His+FcRn(B2M) at 2 μg/mL can bind Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab)

Anti-FcRn (FCGRT & B2M) Reference Antibody (rozanolixizumab)

For Research Use Only. Not for use in diagnostic procedures.