

Anti-Semaphorin 3D (SEMA3D) Mouse Monoclonal Antibody [Clone ID: OTI2H4]

Catalog Number: M09239


Overview

Product Name	Anti-Semaphorin 3D (SEMA3D) Mouse Monoclonal Antibody [Clone ID: OTI2H4]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio SEMA3D mouse monoclonal antibody, clone OTI2H4. Catalog# M09239. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI2H4
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O95025


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 200-291 of human SEMA3D (NP_689967) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IHC 1:500


Anti-Semaphorin 3D (SEMA3D) Mouse Monoclonal Antibody [Clone ID: OTI2H4] (M09239) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SEMA3D (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SEMA3D (1:500).



Immunohistochemical staining of paraffin-embedded Human embryonic brain cortex tissue within the normal limits using anti-SEMA3D mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Human placenta tissue within the normal limits using anti-SEMA3D mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Semaphorin 3D (SEMA3D) Mouse Monoclonal Antibody [Clone ID: OTI2H4]

For Research Use Only. Not for use in diagnostic procedures.