

## Anti-ANGPTL1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]

Catalog Number: M09270


### Overview

Product Name	Anti-ANGPTL1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]
Reactive Species	Human, Mouse
Description	Boster Bio ANGPTL1 mouse monoclonal antibody, clone OTI1B2 (formerly 1B2). Catalog# M09270. Tested in IHC, WB. This antibody reacts with Human, Mouse.
Application	IHC, WB
Clonality	Monoclonal OTI1B2
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O95841


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 24-491 of human ANGPTL1 (NP_004664) produced in E.coli.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:2000 IHC: 1:150


## Anti-ANGPTL1 Mouse Monoclonal Antibody [Clone ID: OTI1B2] (M09270) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ANGPTL1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ANGPTL1.



Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-ANGPTL1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-ANGPTL1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ANGPTL1 Mouse Monoclonal Antibody [Clone ID: OTI1B2]

For Research Use Only. Not for use in diagnostic procedures.