

## Anti-ROR beta Rabbit Monoclonal Antibody

Catalog Number: M09711-1


### Overview

Product Name	Anti-ROR beta Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ROR beta Rabbit Monoclonal Antibody catalog # M09711-1. Tested in WB application. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal 22R66
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92753


### Technical Details

Immunogen	A synthesized peptide derived from human ROR beta
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

## Anti-ROR beta Rabbit Monoclonal Antibody (M09711-1) Images



Western blot analysis of ROR beta expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

### 1 Publications Citing This Product

1. PubMed ID: 26064167, Gong T, Wang Cf, Yuan Jr, Li Y, Gu Jf, Zhao Bj, Zhang L, Jia Xb, Feng L, Liu Sl. Evid Based Complement Alternat Med. 2015;2015:630760. Doi: 10.1155/2015/630760. Epub 2015 May 3. Inhibition Of Tumor Growth And Immunomodulatory Effects Of Flavonoids...

Visit [bosterbio.com/anti-ror-beta-rabbit-monoclonal-antibody-m09711-1-boster.html](http://bosterbio.com/anti-ror-beta-rabbit-monoclonal-antibody-m09711-1-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ROR beta Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.