

Anti-TPIPb Antibody (C-term)(Ascites)

Catalog Number: M09798


Overview

Product Name	Anti-TPIPb Antibody (C-term)(Ascites)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TPIPb Antibody (C-term)(Ascites) (Catalog # M09798). Tested in WB application(s). This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal 978CT6.1.4
Formulation	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q6XPS3

Technical Details

Immunogen	Purified His-tagged TPIPb protein was used to produced this monoclonal antibody.
Isotype	IgG2a
Suggested Dilutions	WB: 1:5000

Anti-TPIPb Antibody (C-term)(Ascites) (M09798) Images



TPIPb Antibody (C-term) western blot analysis in 293,A431,A2058,mouse NIH/3T3 cell line lysates (35ug/lane).This demonstrates the TPIPb antibody detected the TPIPb protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TPIPb Antibody (C-term)(Ascites)

For Research Use Only. Not for use in diagnostic procedures.