

Anti-Pepsin A-5 PGA5 Monoclonal Antibody

Catalog Number: M10094


Overview

Product Name	Anti-Pepsin A-5 PGA5 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Pepsin A-5 PGA5 Monoclonal Antibody catalog # M10094. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 2C1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P0DJ9


Technical Details

Immunogen	Purified recombinant fragment of human PGA5 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-Pepsin A-5 PGA5 Monoclonal Antibody (M10094) Images



Western blotting validation for Anti-Pepsin A-5 PGA5 Monoclonal Antibody M10094 Western Blot (WB) analysis using PGA5 Monoclonal Antibody against HepG2 (1) and SMMC-7721 (2) cell lysate. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of PGA5 using Anti-Pepsin A-5 PGA5 Monoclonal Antibody (M10094). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Stomach cancer tissues with DAB staining using PGA5 Monoclonal Antibody. For more protocol information of IHC

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pepsin A-5 PGA5 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.